

minicap sebia flex piercing



Small footprint
for full features



©2020 **sebia**  Non-contractual illustration - 03/2018 - Ref.: MFPPB1EN - Shutterstock: ©Borislav Borisov
Sebia reserves the right to modify at any time and without notice the information contained in this document intended for healthcare professionals.
This document contains information on products which is targeted to a wide range of audiences and could contain product details or information otherwise not accessible or valid in your country.


Benelux - Tel.: +32 2 702 64 64 (BE) - Tel.: +31 20 610 8099 (NL) - E-mail: sebia.benelux@sebia.com
Brasil - Tel.: +55 11 3849 0148 - E-mail: sebia@sebia.com.br
China - Tel.: +86 (21) 6360 1436 - Fax: +86 (21) 6361 2011 - E-mail: china@sebia.com
Dubai - Tel.: +971 (0)4 878 6977
Germany - Tel.: +49 (0) 661 33081 - Fax: +49 (0) 661 31881 - E-mail: sebia@sebia.de
India - Tel.: +91 124 4262100 - E-mail: sebiaindia@sebia.com
Italy - Tel.: +39 055 24851 - Fax: +39 055 0982083 - E-mail: info@sebia.it
Portugal - Tel.: +351 211 224 653 - Fax: +351 211 305 079
UK - Tel.: +44 (01276) 600636 - Fax: +44 (01276) 38827 - E-mail: sales@sebia.co.uk
Spain - Tel.: +34 93 208 15 52 - Fax: +34 93 458 55 86 - E-mail: sebia@sebia.es
Swiss - Tel.: +41 44 787 88 10 - Fax: +41 44 787 88 19 - E-mail: contact.ch@sebia.com
U.S.A. - Tel.: +1 (770) 446-3707 - Fax: +1 (770) 446-8511 - E-mail: info@sebia-usa.com

sebia

RCS Evry 672 041 902
Parc technologique Léonard de Vinci
27 rue Léonard de Vinci
CP 8010 Lisses - 91008 Evry cedex
France
Tel. : +33 1 69 89 80 80
E-mail : sebia@sebia.com
www.sebia.com

sebia

minicap sebia flex piercing



MINICAP Flex Piercing

Multi-parameter instrument, Large test menu

- Serum & Urine Proteins Electrophoresis (SPE & UPE), Serum & Urine Immunotyping, Hb A1c (venous blood) Hemoglobin electrophoresis (adult whole blood), CDT/CDT_{IFCC}

Analytical performances

- Based on a solid experience of more than 10,000 capillary electrophoresis systems worldwide

Full automation

- Sample mixing and cap piercing capabilities on Hb A1c & Hemoglobin technics. Analytical superiority and enhanced results precision. Maximizes productivity and improve operator's safety.






- Positive sample identification:** full traceability from primary tube to final result

Throughput

- 20 tests/h (SPE & UPE)
- 8 tests/h (Hb A1c & Hemoglobin)
- 10 tests/h (CDT & CDT_{IFCC})

Specifications

- Small footprint
- Continuous sample loading




Main functions of the SQL Client/Server database

- Networking of all Sebia instruments
 - within a laboratory
 - between off-site locations
- Remote validation of the results
- Unlimited data storage
- Web access


Interpretation software and results validation

- Intuitive and user-friendly interface
- Multicriteria search function
- Automatic recall of patient's history
- Automatic triage of the results



Technology advancements

Sample mixing



Cap piercing

